# NOOA

## BODY ANALYSER SCALE

### Instruction Manual

Model:NO601H

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Thank You for purchasing your new NOAA Body Analyser Scale.

We recommend that you spend some time reading this instruction manual in order to fully understand all the operational features it offers. This scale is designed to analyse proportion of body fat, body water, body muscles and bone content. With reference to your gender, age, height and weight, the scale is able to determine your daily calorie requirement for reference.

Read all the safety instructions carefully before use and keep this instruction manual for future reference.

## ■ GENERAL TIPS

### 1. During Measuring

- A. The device is designed to allow auto step on. Please always ignore the first reading and only start taking effective record from the second time weighing.
- B. Always remove your shoes and socks, and clean your feet before use.
- C. Step on the weighing platform gently.
- D. To attain a compatible body analysis, it is suggest to weigh yourself at the same time of day (preferably in the morning) after visiting the toilet, with an empty stomach and without clothing.
- E. The analysis may not be accurate if you weigh yourself after intensive exercise, intensive diet or at extreme dehydration condition.
- F. Place the scales on a solid, flat horizontal surface. Do not place on a carpet or any soft surface. Do not move during weighing.
- G. The body fat /water content analysed by the scale may have deviation for the person who are:
  - children aged under 10 or adults aged over 100 ( They can only use the normal weighing mode)
  - body building athletes or other occupational athletes.
- H. The product is not suitable for use by person who is pregnant, with symptom of edema, with a heart pacemaker or other implanted medical device, or on dialysis treatments.

### 2. Safety Warning

- A. The body fat, water content and muscle content analysed by the scale is for reference only. If your body analysis result is over/ below the normal level, always consult with your doctor for more advice.
- B. This device has delicate electronic components. Therefore interference may occur when it is near devices which transmit radio signals. If errors occur in the display, remove away from these signals or remove the battery briefly and then insert it again.
- C. Do not step onto the scale with wet feet and do not use the scale if the platform is moist. Otherwise you could slip.

## ■ FEATURES / SPECIFICATIONS

1. LED body fat level indicator
2. Using high precision strain gauge sensor
3. Analysing proportion of body fat, body water, body muscles and bone and calorie intake requirement
4. Built-in memory for 12 different users
5. Auto step on weighing and consecutive add-on weighing without entering personal data
6. Low battery and overload indication
7. Memory recall from last user
8. Compare current result to the last record
9. Analyse body fat condition into 5 levels


## ■ PERFORMANCE INDEX

Item Index	Weighing Capacity	Fat%	Hydration%	Muscle%	Bone%	Calorie
Index Range		4.0% ~60.0%	27.5% ~66.0%	20.0% ~56.0%	2.0% ~20.0%	
Division	0.1 kg 0.2 lb	0.1 %	0.1 %	0.1 %	0.1 %	1KCAL

## ■ REPLACING BATTERIES


Depending on the model of your scale, different battery type is required. Find the battery type of your scale, and replace the battery as follows:

- (A) The scale is operated by one CR2032 (3V) lithium battery (included). Before first use, remove the insulation strip at the battery compartment.
- (B) The scale is operated by two AA or AAA (1.5V) batteries (not included).

PIC 1

1. Open the battery compartment cover.
2. Remove the old batteries (if any) using a sharp object as shown in the picture. Install one CR2032 (3V) lithium battery. Make sure the positive (+) is facing upwards. HINT: Make sure one side of the battery is placed underneath the battery baffle, then push the other side of the battery until it clicks and secured.
3. Close the battery compartment cover.

PIC 2

1. Open the battery compartment cover.
2. Remove the old batteries (if any). Insert two AA or AAA (1.5V) batteries (not included). Make sure to match the (+) and (-) ends of the batteries with the (+) and (-) ends indicated in the battery compartment.
3. Close the battery compartment cover.

## FUNCTION BUTTONS


**SET ON/(SET) :** Press to switch the scale on.

Press to confirm the value entered in setting personal profile.

- ▲ : Press the Up/ Down button to select your desired value in setting personal profile.
- ▼ : Press the Up/ Down button to select your desired value in setting personal profile.


**HINT:** Press and hold the Up/ Down button to speed up the process of setting the values.

## CHANGING WEIGHT UNIT

Depending on the model of your scale, your scale may support weight units in kg/lb/st.


Press the UNIT button repeatedly at the rear of the scale to select your desired unit of weight. The units will be displayed in the LC display.

After the weight unit is confirmed, the scale will automatically switch off in 3 seconds. The selected unit will be shown on the display when the scale is switched on again.

**HINT:** If the weight unit "lb/st" is selected, "feet/inch" is the default height unit; if the weight unit "kg" is selected, "cm" is the default height unit.

## LC DISPLAY

## ACTIVITY LEVEL

Inactive

Light activity

Average activity


Above average activity

Very active

## NORMAL WEIGHING MODE (without entering personal data)

### (1) Auto step on weighing


1. Place the scales on a solid, flat horizontal surface. Do not place on a carpet or any soft surface. Always make sure there are no other objects beneath the scale before weighing.
2. Place even weight on each leg and do not move. Moving during the weighing process can lead to discrepancies in the results. The scales will start to measure immediately. After flashing, the weight display will show the result.
3. Step off the platform and view the weight measurement.

## (2) Consecutive add-on weighing

Consecutive add-on weight can also be measured after the first weight shown on the display.

NOTE: The add-on weight should exceed 2kg.

HINT: This function is useful if you want to weigh your baby. See below procedures for details.

1. First weigh yourself by placing even weight on each leg and do not move. The scales will start to measure immediately. After flashing, your weight will be shown on the display.
2. Do not step off the platform. Now hold your baby and weigh together. The scales will start to measure immediately. After flashing, the new weight measurement will be displayed.
3. Step off the platform and view the new weight measurement. By subtracting your own weight from the new weight measurement, you will get the weight of your baby.

## PERSONAL BODY ANALYSIS MODE

### (1) Entering personal data

① First time use, the personal data settings are defaulted as below:

Gender	Height	Activity level	Age	Target weight
Male	165cm/5'05.0"		25	50kg/110lb


② Personal data's range can be set as below:

User ID code	Gender	Height	Activity level	Age	Target weight
1~12	male/female	100~250cm 3'03.5"~8'02.5"		10~100	20~Max. Capacity

### ③ Entering personal data

1. Press the SET (ON/SET) button to switch the scale on and to set your personal data.  
Note: The User ID will flash in the LC display. The number stands for the data profile of one person from the user group.
2. Press the Up/ Down button to select the data profile (1–12). Press the SET (ON/SET) button to confirm the selected data profile.
3. Press the Up/ Down button to select the option of the gender. Press the SET (ON/SET) button to confirm the selected data profile.
4. Press the Up/ Down button to select your activity level. Press the SET (ON/SET) button to confirm the selected data profile.
5. Press the Up/ Down button to select your height. Press the SET (ON/SET) button to confirm the selected data profile.
6. Press the Up/ Down button to select your age. Press the SET (ON/SET) button to confirm the selected data profile.
7. Press the Up/ Down button to select your target weight.  
After you have made all these settings, gently tap on the scale platform with your hand or foot.
8. The display will show zero and the scale will switch off automatically after about 8 seconds if it is not used. Your personal data have successfully been saved. You can either make settings for further users or measure your body weight.

NOTE: If the display show zero during the personal profile setting, you can always press the SET (ON/SET) button to restart the setting mode and re-enter the User ID and personal data again.  
HINT: When setting the target weight, you may press and hold the Up/ Down button to speed up the process of setting the values with increment of 1kg.


## (2) Measuring the body weight and analysis of proportion of body fat, body water and body muscles.


Always remove your shoes and socks, and clean your feet before use. Place even weight on each leg and do not move. Make sure your feet are in good contact with the two electrodes on the scale platform.

1. Press the SET (ON/SET) button, the display will show the last measurement record of the last User ID code.
2. Press the Up/ Down button to select the data profile (1–12). Gently tap on the scale platform with your hand or foot. The display will show the last measurement record of your User ID code.
3. Wait until the display show zero. Step barefoot onto the platform and do not move. Ensure that you are standing on the two electrodes. The scales will start to measure immediately. After flashing, the weight display will show the result.
4. The body fat/ body water/ body muscle/ bone/ calorie analysis is beginning; remain standing on the two electrodes. The measurement lasts a few seconds (display indication oooo) and the body fat index LED bar at the middle of the scale will light on.
5. The display then shows the measured values three times to show current weight, increased / decreased weight compared to last measurement, body fat index, the percentage of body fat , body water, body muscle, bone and calorie requirement. The scales switch off automatically after these values have been displayed.

NOTE: The display will show "C", if you do not stand on the two electrodes during weighing.

## WARNING INDICATIONS

### Low Battery warning

If the scale shows "Lo", it indicates that the batteries are too weak and it is time to replace the batteries. Refer to the section "REPLACING BATTERIES" to replace batteries.


### Overload warning

If the scale shows "Err", it indicates that the scale is overloaded. The maximum capacity is indicated on the scale. Step off the scale to avoid causing damage to the scale.


### Low Fat warning

If the scale shows "ErrL", it indicates that your body fat portion is too low. Please consult your doctor or nutritionist for professional advice.

### High Fat warning

If the scale shows "ErrH", it indicates that your body fat portion is too high. Please consult your doctor or nutritionist for professional advice.

### Error warning

If the scale shows "C", it indicates that there are measuring errors. Please step off the scale and retry.

## ■ ABOUT FAT

### 1. About Fat%

Fat % is a % index of fat content inside a human body. The fat content in a human body is in a stable level. Any big fluctuations in fat will have a bad effects on the metabolism, which may even cause various illnesses. Hence it is always very important to keep track on our fat% in order to have a healthy body.

### 2. Working Principles

By applying a low-voltage electrical current of a safe frequency through the human body, with the fact that the impedance of fat is bigger than that of muscle or other body issues, this body analyser scale can analyse the biological resistance which reflect the fat% in the body.

### 3. Fat/Hydration Content Fitness Assessment Chart

After attaining the analysed result, you may use the below assessment table for reference.

Age	Fat(F)	Hydration(F)	Fat(M)	Hydration(M)	Status	Body fat index
≤30	4.0-16.0%	66.0-57.8%	4.0-11.0%	66.0-61.2%	Underweight	
	16.1-20.5%	57.7-54.7%	11.1-15.5%	61.1-58.1%	Thin	
	20.6-25.0%	54.6-51.6%	15.6-20.0%	58.0-55.0%	Normal	
	25.1-30.5%	51.5-47.8%	20.1-24.5%	54.9-51.9%	Fat	
	30.6-60.0%	47.7-27.5%	24.6-60.0%	51.8-27.5%	Overweight	
>30	4.0-20.0%	66.0-55.0%	4.0-15.0%	66.0-58.4%	Underweight	
	20.1-25.0%	54.9-51.6%	15.1-19.5%	58.3-55.3%	Thin	
	25.1-30.0%	51.5-48.1%	19.6-24.0%	55.2-52.3%	Normal	
	30.1-35.0%	48.0-44.7%	24.1-28.5%	52.2-49.2%	Fat	
	35.1-60.0%	44.6-27.5%	28.6-60.0%	49.1-27.5%	Overweight	

• The above information is for reference only.

## ■ ABOUT CALORIE

### 1. Calorie & Energy Consumption

Calorie is a unit for energy. The calorie requirement estimated by this scale indicates the amount of energy that you required to sustain a normal daily metabolism and activities with reference to your weight, height, age and gender.

### 2. Calorie & Losing Weight

Exercising is an effective way to lose weight by controlling calorie intake. When the daily energy intake cannot support the energy consumption of the human body, the fat which are stored inside the body will decompose, in order to provide energy to the human body, as a result this lead to weight loss.

## ⚠ CARE AND MAINTENENCE

1. Do not step onto the scale with wet feet and do not use the scale if the platform is moist.
2. Keep standing still during the course of weighing.
3. DO NOT strike, shock or drop the scale.
4. Clean the scale using a slightly moist cloth. Do not use chemical cleaners or abrasive cleaning agents for cleaning the scale. Ensure there is no water entering the scale. If this does occur, then wait until the device has dried out fully before you use it again.
5. Keep the scale in a cool dry area.
6. Always keep you scale in a horizontal position.
7. If the scale does not appear to work, check if the batteries have been inserted correctly and reinsert them if necessary.  
If the scale shows “Lo”, replace with new batteries.
8. This device has delicate electronic components. Therefore interference may occur when it is near devices which transmit radio signals. If errors occur in the display, remove away from these signals or remove the battery briefly and then insert it again.
9. The scale is intended only for use in a private domestic environment and not for commercial use.

## ■ ABOUT FAT

### 1. About Fat%

Fat % is a % index of fat content inside a human body. The fat content in a human body is in a stable level. Any big fluctuations in fat will have a bad effects on the metabolism, which may even cause various illnesses. Hence it is always very important to keep track on our fat% in order to have a healthy body.

### 2. Working Principles

By applying a low-voltage electrical current of a safe frequency through the human body, with the fact that the impedance of fat is bigger than that of muscle or other body issues, this body analyser scale can analyse the biological resistance which reflect the fat% in the body.

### 3. Fat/Hydration Content Fitness Assessment Chart

After attaining the analysed result, you may use the below assessment table for reference.

Age	Fat(F)	Hydration(F)	Fat(M)	Hydration(M)	Status	Body fat index
≤30	4.0-16.0%	66.0-57.8%	4.0-11.0%	66.0-61.2%	Underweight	
	16.1-20.5%	57.7-54.7%	11.1-15.5%	61.1-58.1%	Thin	
	20.6-25.0%	54.6-51.6%	15.6-20.0%	58.0-55.0%	Normal	
	25.1-30.5%	51.5-47.8%	20.1-24.5%	54.9-51.9%	Fat	
	30.6-60.0%	47.7-27.5%	24.6-60.0%	51.8-27.5%	Overweight	
>30	4.0-20.0%	66.0-55.0%	4.0-15.0%	66.0-58.4%	Underweight	
	20.1-25.0%	54.9-51.6%	15.1-19.5%	58.3-55.3%	Thin	
	25.1-30.0%	51.5-48.1%	19.6-24.0%	55.2-52.3%	Normal	
	30.1-35.0%	48.0-44.7%	24.1-28.5%	52.2-49.2%	Fat	
	35.1-60.0%	44.6-27.5%	28.6-60.0%	49.1-27.5%	Overweight	

• The above information is for reference only.

## ■ ABOUT CALORIE

### 1. Calorie & Energy Consumption

Calorie is a unit for energy. The calorie requirement estimated by this scale indicates the amount of energy that you required to sustain a normal daily metabolism and activities with reference to your weight, height, age and gender.

### 2. Calorie & Losing Weight

Exercising is an effective way to lose weight by controlling calorie intake. When the daily energy intake cannot support the energy consumption of the human body, the fat which are stored inside the body will decompose, in order to provide energy to the human body, as a result this lead to weight loss.

## ⚠ CARE AND MAINTENENCE

1. Do not step onto the scale with wet feet and do not use the scale if the platform is moist.
2. Keep standing still during the course of weighing.
3. DO NOT strike, shock or drop the scale.
4. Clean the scale using a slightly moist cloth. Do not use chemical cleaners or abrasive cleaning agents for cleaning the scale. Ensure there is no water entering the scale. If this does occur, then wait until the device has dried out fully before you use it again.
5. Keep the scale in a cool dry area.
6. Always keep you scale in a horizontal position.
7. If the scale does not appear to work, check if the batteries have been inserted correctly and reinsert them if necessary.  
If the scale shows “Lo”, replace with new batteries.
8. This device has delicate electronic components. Therefore interference may occur when it is near devices which transmit radio signals. If errors occur in the display, remove away from these signals or remove the battery briefly and then insert it again.
9. The scale is intended only for use in a private domestic environment and not for commercial use.

This symbol on the product or in the instructions means that your electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There are separate collection systems for recycling in the EU.

For more information, please contact the local authority or your retailer where you purchased the product.

Imported and Exclusively marketed by:  
Expert AS  
PO Box 43,  
N-1483 Hagan

Customer helpline number Norge: 815 48 100  
Customer helpline number Finland: 020 7100 670  
Customer helpline number Danmark +45 70 70 17 07  
Email Norge: [kundesenter@expert.no](mailto:kundesenter@expert.no)  
Email Finland: [asiakaspalvelu@expert.fi](mailto:asiakaspalvelu@expert.fi)  
Email Danmark: [kundeservice@expert.dk](mailto:kundeservice@expert.dk)

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