NOOA

PRO PERFORMANCE

Body Analyser Scale

Download on the
App Store

Requires iOS 8.0 or later

ANDROID APP ON
Google play

Requires Android 4.3 or later

NOBS90BT

User Manual

GB

Thank you for purchasing your new NOOA body analyser scale. These operating instructions will help you use it properly and safely. We recommend that you spend some time reading this instruction manual in order that you fully understand all the operational features it offers. Read all the safety instructions carefully before use and keep this instruction manual for future reference.

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SAFETY WARNING

1. Read all the instructions carefully before using this attachment and keep them for future reference.
2. This device has delicate electronic components. Therefore interference may occur when it is near devices which transmit radio signals. If errors occur in the display, remove away from these signals or remove the battery briefly and then insert it again.
3. Keep batteries away from children. Do not leave batteries lying around. They could be swallowed by children or pets. In the event of swallowing a battery, contact a doctor immediately.
4. This scale is for household use only. Do not use for commercial purpose.
5. The glass platform will be slippery when it is wet. Keep it dry!
6. Keep standing still during the course of weighing.
7. Do not strike, shock or drop the scale.
8. Treat the scale with care as it is a precise instrument.
9. Clean the scale with a damp cloth and avoid water getting inside the scale. Do not use chemical/abrasive cleaning agents.

10. Keep the scale in cool place with dry air.
11. Always keep your scale in horizontal position.
12. If failed to turn on the scale, check if battery is installed, or battery power is low. Insert or replace new battery.
13. If there is error on display or can not turn off the scale for a long time, please take off the battery for about 3 seconds, then install it again to remove soft-fault. If you can not resolve the problem, please contact your dealer for consultation.

PRODUCT OVERVIEW

PERFORMANCE INDEX

| Item Index | Weighing Capacity | Fat% | Hydration% | Muscle% | Bone% | Calorie |
|-------------|-------------------|----------------|-----------------|-----------------|----------------|---------|
| Index Range | | 4.0% ~60.0% | 27.5% ~66.0% | 20.0% ~56.0% | 2.0% ~20.0% | |
| Division | 0.1kg | 0.1% | 0.1% | 0.1% | 0.1% | 1KCAL |

BATTERY INSTALLATION


| Name | Illustration | Installation |
|---------------------------------------|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2x1.5V AAA Battery | | AAA battery installation Place install battery according to polarity. |
| 1x3V or 2x3V Lithium battery (CR2032) | | A. Take out the used battery with the aid of a cuspidal object. B. Install new battery by putting one side of the battery beneath the battery baffle and then pressing down the other side. Please take note to the polarity. |

Note: Please use corresponding battery. Please remove battery if it is not in use in a long time (more than three months).

KEY ILLUSTRATION

- SET** : Turn on the scale and prompt to the next parameter.
(Please refer to the Setting Personal Parameter overleaf)
- ▲** : Value up by each interval press.
Press and hold the button to have consecutive value added.
- ▼** : Value down by each interval press.
Press and hold the button to have consecutive value reduced


LCD DISPLAY


TO CONVERT WEIGHT UNIT

If the scale is with both measuring system of Metric (kilograms) and Imperial (pounds, stones), you can convert the weight unit as follows:

- Condition 1: if the unit switch is on the back of the scale,
- Find the unit switch on the back of scale.
 - Pull up and down the unit switch to convert the scale unit between kg-lb-st.


- Condition 2: if the unit switch is on the back of the scale,
- Find the unit switch on the back of the scale.
 - Press the unit conversion button when the scale is powered on, and the LCD shows current weight unit. Press the button again to convert the weight unit.

NORMAL WEIGHING MODE (by passing parameter setting)

1. Auto step on weighing

After installing the battery, place the scale on a flat floor. Stand still on scale when testing. Then wait for seconds, display will show the same weight reading for approximately a few seconds.


Install the battery then place the scale on a flat hard surface to ensure an accurate reading.

Lightly step on the scale (Auto-On) then wait for a reading.

The reading result will flash after the weight has stabilised.

2. Consecutive add-on weighing

After the digits have stabilized, additional loadings can be conducted. Ensure each add-on loading exceeds 2kg to get a new reading. For example: If you want to get a weight value of your baby, you can weigh yourself first, and then hold your baby together, without stepping off the scale. Then you will get the total weight reading. Subtracting the previous reading of yourself from the new reading, you will get the difference (meaning the weight of your baby)..

Digits flash, stabilized and locked

Display total weight.

FAT/HYDRATION/MUSCLE/BONE/CALORIE MEASURING MODE

(1) Setting personal parameters

A. The first time you install the battery and switch on, all default parameter will be set as follows:


| Memory | Gender | Height | Age |
|--------|--------|---------------|-----|
| 1 | Female | 165cm/5'05.0" | 25 |

B. Parameters can be set in the following range:

| Memory | Gender | Height | Age |
|--------|-------------|------------------------------|--------|
| 1~12 | male/female | 100~250cm 3'03.5"-8'02.5" | 10~100 |

C. Parameter setting

Example: the 3rd user, female, 168cm, 28 years old, can set her parameters as follow:


NOTE:

- During parameter setting, if idle for 6 seconds, the parameters already set will be memorized automatically and the LCD shows "0", meaning that body fat measuring will be started.
- If the LCD shows "0" during the parameter setting, press SET for 3 seconds to return to the setting mode, then you can set the user code, gender, height and age

(2) To Start Measuring (Make sure your feet are touching with the electrodes on the scale platform, or your body fat can not be measured.)

To ensure accuracy, please remove your shoes & socks before stepping on.
Example: the 3rd user, female, 168cm, 28 years old can measure as follow:

WARNING INDICATIONS

Low Battery Indication

The battery power is running low, please replace with a new battery.

Over-load Indication

The weighing subject on the platform exceeds the maximum capability of scale. Please step off to avoid damage.

Low Fat% Indication

The Fat% is too low. Please lead a more nutritional diet and take good care of yourself.

High Fat% Indication

The Fat% is too high. Please watch your diet and do more exercise.

Re-test

Error exists, please re-test to get the correct result.

ABOUT FAT

1. About Fat%

- Fat% is an index % of fat content in human body.
- The fat content in the human body is of a certain level, overabundance of fat will bring bad effects on metabolism, thereby, various illnesses may occur.

2. Working Principles

As electrical signals of a certain frequency transmit through the human body, the impedance of fat is stronger than that of muscle or other human tissues. Through the use of biological impedance analysis, this product will pass a safe frequency of low voltage electrical current through human body, thus, the Fat% will be reflected according to the resistance differences when the electrical signals are transmitting through human body.

3. Fat/Hydration Content Fitness Assessment Chart

After weighing yourself, you can use the following Assessment Chart for reference:

| Age | Fat(F) | Hydration(F) | Fat(M) | Hydration(M) | Status | Key Line Indication |
|-----|------------|--------------|------------|--------------|------------|---------------------|
| ≤30 | 4.0-16.0% | 66.0-57.8% | 4.0-11.0% | 66.0-61.2% | Thin | |
| | 16.1-20.5% | 57.7-54.7% | 11.1-15.5% | 61.1-58.1% | Slim | |
| | 20.6-25.0% | 54.6-51.6% | 15.6-20.0% | 58.0-55.0% | Normal | |
| | 25.1-30.5% | 51.5-47.8% | 20.1-24.5% | 54.9-51.9% | Heavy | |
| | 30.6-60.0% | 47.7-27.5% | 24.6-60.0% | 51.8-27.5% | Overweight | |
| >30 | 4.0-20.0% | 66.0-55.0% | 4.0-15.0% | 66.0-58.4% | Thin | |
| | 20.1-25.0% | 54.9-51.6% | 15.1-19.5% | 58.3-55.3% | Slim | |
| | 25.1-30.0% | 51.5-48.1% | 19.6-24.0% | 55.2-52.3% | Normal | |
| | 30.1-35.0% | 48.0-44.7% | 24.1-28.5% | 52.2-49.2% | Heavy | |
| | 35.1-60.0% | 44.6-27.5% | 28.6-60.0% | 49.1-27.5% | Overweight | |

The above information is only for reference

ABOUT CALORIE

1. Calorie & Energy Consumption

Calorie: Energy Unit

The data reading indicates the energy you need to sustain the daily metabolism and ordinary activities required for your weight, height, age and sex.

2. Calorie & Losing Weight

Controlling your calorie intake combined with proper exercise is an effective way to lose weight. This is a very simple theory. As the daily energy intake can not support the energy consumption of the human body, the saccharide and fat stored inside will decompose and provide energy to the human body. This will leads to weight loss.

APP INSTRUCTION

- How to download the application
 - Scan the QR code on the bottom of scale to download the Application.
 - At the "App Store" or "Google Play", search the name of App "BodyMonitor" to download the application.

- How to connect Bluetooth
 - Make sure that Bluetooth is turned on before use
 - If not enabled, go to the "Settings Menu" to turn it on.
 - Press the scale or stand on it.

- How to create a new user
 - Click the app icon, to enter "user settings".
 - Complete the information to enter weighing interface.

- How to use the scale with the app
 - Press scale to turn it on
 - Step on the scale.
 - Your weight will be shown on the scale and the App. After several seconds your body fat, hydration, muscles, bone percentage and kcal suggestions will be shown.
 - All measurements will be saved automatically to data history.
 - You can review the history, add new users and set reminders for your weighing habit.

- Requirements for mobile devices:
 - Bluetooth version: Bluetooth 4.0 or above
 - Apple System: iPhone 4S or later, iOS 8.0 or later
 - Android System: Android 4.3 or later

Trouble Shooting

If the software can't obtain weighing results from Bluetooth scale, possible solutions could be as blow:

- Confirm the connection sign on the App interface to know the connection status.
- Confirm Bluetooth is on in your device.
- Confirm Bluetooth scale is on within the valid range of signal coverage.
- Relaunch the application.
- Turn off Bluetooth and on again, then bridge the Bluetooth connection between scale and smart phone again.

Failure to connect to your device

Please make sure only one device is set up for connection with the Bluetooth scale. Please make sure your device is equipped with Bluetooth 4.0 or above.

LCD has no response

First check that the batteries are still live, if not, replace them if necessary.
Repairs may only be carried out by our authorized customer service or dealers.

ADVISE FOR USE & CARE

1. DO NOT use the scale in humid, overheat or overcold environment. (Temperature range form 5°C to 35°C)
2. DO NOT use the scale on soft mat, lawn, slippery floor or uneven floor.
3. DO NOT drop, shock or strike the scale as it is a precise instrument.
4. Keep the scale clean. DO NOT use acid cleaning agent.
5. DO NOT use the scale on a place with intense shake. Keep standing firmly on scale when in use.
6. To ensure accuracy, please test at the same time of the sequential days.

FEATURES/SPECIFICATIONS

A. Scale Body

1. Model: NOBS90BT
2. 74 x 35,6 mm LCD
3. Low battery / over load indication
4. Capacity: 180 kg / 400 lb / 28st
5. Graduation: 100g / 0.2 lb
6. Power: 2 x 1.5V AAA
7. Uses a high precision strain gauge sensor
8. Multifunction: tests body fat %, body hydration %, body muscle%, body bone% and suggests calorie intake requirements.
9. Built-in memory for 12 different users
10. Optional normal weighing mode (bypassing parameter setting): Auto step on weighing and consecutive add-on weighings
11. Low battery and overload indication

B. Bluetooth

1. Directly check personal weight, body fat%, body hydration%, body muscle%, and body bone% information through graphs.
2. Set password to protect personal information and store multi-user data;
3. Simple and clear graphs for current weight trends, form which you are able to keep track of changes;
4. Bluetooth wireless technology;
5. Connect with different smart devices respectively;

This symbol on the product or in the instructions means that your electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There are separate collection systems for recycling in the EU.

Notes for battery disposal

The batteries used with this product contain chemicals that are harmful to the environment. To preserve our environment, dispose of used batteries according to your local laws or regulations. Do not dispose of batteries with normal household waste.

For more information, please contact the local authority or your retailer where you purchased the product.

Imported and Exclusively marketed by:
Power International AS,
PO Box 523, N-1471 Lørenskog, Norway

Power Norge:
<https://www.power.no/kundeservice/> T: 21 00 40 00
Expert Danmark:
<https://www.expert.dk/kundeservice/> T: 70 70 17 07
Power Danmark:
<https://www.power.dk/kundeservice/> T: 70 33 80 80
Punkt 1 Danmark:
<https://www.punkt1.dk/kundeservice/> T: 70 70 17 07
Expert Finland:
<https://www.expert.fi/tuki/asiakaspalvelu/> T: 020 7100 670
Power Finland:
<https://www.power.fi/tuki/asiakaspalvelu/> T: 0305 0305
Power Sverige:
<https://www.power.se/kundservice/> T: 08 517 66 000

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