NOOA

PRO PERFORMANCE

Body Analyser Scale

NOPS972
User Manual

GB

Thank you for purchasing your new SENZ body analyser scale. These operating instructions will help you use it properly and safely. We recommend that you spend some time reading this instruction manual in order that you fully understand all the operational features it offers. Read all the safety instructions carefully before use and keep this instruction manual for future reference.

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SE

Tack för att du köpt en ny SENZ kroppsanalysvåg. These operating instructions will help you use it properly and safely. We recommend that you spend some time reading this instruction manual in order that you fully understand all the operational features it offers. Read all the safety instructions carefully before use and keep this instruction manual for future reference.


| | | |
|------------------------------------|------|---------|
| Instruction Manual - English | page | 4 - 15 |
| Bruksanvisning- Norsk | sida | 16 - 27 |
| Käyttöohjeet - Suomi | side | 28 - 39 |
| Brugsanvisning - Dansk | side | 40 - 51 |
| Bruksanvisning - Svenska | sivu | 52 - 63 |

SAFETY WARNING

1. Always weigh / use the scale on a solid and flat surface. The device is designed to allow auto step on. Always ignore the first reading and only start taking effective records from the second weighing.
2. Always remove your shoes and socks , and wipe clean your feet before using.
3. Step on the weighing platform gently.
4. It is advisable to take your measurements at the same time of the day.
5. Measured result could be misleading after intensive exercise, excessive dieting or under extreme dehydration condition.
6. Keep away from strong electro-magnetic field when using the scale.
7. Children under 10 or adults over 100 can only use the device under normal weighing mode.
8. The measured data of the following people may have some deviation:
 - Adults over 70
 - Body building athletes or other occupational athletes.
9. People with the following conditions should not use the device:
 - Pregnant women, Those with symptoms of edema, or on dialysis treatment.

10. The Fat% Hydration% content evaluated by the Fat Monitor Scale is only for reference (not for medical purpose). If your Fat% Hydration% content is over/below normal level, please consult your doctor for more advice.
11. If the scale fails to turn on, check whether the battery is installed, or if the battery power is low. Insert or replace new battery.
12. Do not use for commercial purpose.
13. It could be slippery if you weigh your self on a wet surface. If the weighing platform surface is glazed, never step on the edge of the scale and always keep your balance during your weighing.
14. People with heart pacemakers or other implanted medical devices should not use the body fat testing function.

PRODUCT OVERVIEW


PERFORMANCE INDEX

| Item Index | Weighing Capacity | Fat% | Hydration% | Muscle% | Bone% | Calorie |
|-------------|-------------------|----------------|-----------------|-----------------|----------------|---------|
| Index Range | | 4.0% ~60.0% | 27.5% ~66.0% | 20.0% ~56.0% | 2.0% ~20.0% | |
| Division | 0.1kg | 0.1% | 0.1% | 0.1% | 0.1% | 1KCAL |

BATTERY INSTALLATION


- Press battery gently, the battery will pop up automatically.
- Install new battery by connecting one side of battery to the positive plate first, and then pressing down the other side. Please observe the polarity.

KEY ILLUSTRATION


- SET** : Turn on the scale and prompt to the next parameter.
(Please refer to the Setting Personal Parameter overleaf)
- ▲ : Value up by each interval press.
Press and hold the button to have consecutive value added.
- ▼ : Value down by each interval press.
Press and hold the button to have consecutive value reduced

LCD DISPLAY

TO CONVERT WEIGHT UNIT


- Find the unit switch on the back of the scale.
- Press the unit conversion button when the scale is powered on, and the LCD shows current weight unit. Press the button again to convert the weight unit.

NORMAL WEIGHING MODE (by passing parameter setting)

1. Auto step on weighing

After installing the battery, place the scale on a flat floor. Stand still on scale when testing. Then wait for seconds, display will show the same weight reading for approximately a few seconds.

Install the battery then place the scale on a flat hard surface to ensure an accurate reading.

Lightly step on the scale (Auto-On) then wait for a reading.

The reading result will flash after the weight has stabilised.

2. Consecutive add-on weighing

After the digits have stabilized, additional loadings can be conducted. Ensure each add-on loading exceeds 2kg to get a new reading.

For example: If you want to get a weight value of your baby, you can weigh yourself first, and then hold your baby together, without stepping off the scale. Then you will get the total weight reading. Subtracting the previous reading of yourself from the new reading, you will get the difference (meaning the weight of your baby)..

FAT/HYDRATION/MUSCLE/BONE/CALORIE MEASURING MODE

(1) Setting personal parameters

A. The first time you install the battery and switch on, all default parameter will be set as follows:


| Memory | Gender | Height | Age |
|--------|--------|---------------|-----|
| 1 | Female | 165cm/5'05.0" | 25 |

B. Parameters can be set in the following range:

| Memory | Gender | Height | Age |
|--------|-------------|------------------------------|--------|
| 1~12 | male/female | 100~250cm 3'03.5"~8'02.5" | 10~100 |

C. Parameter setting

Example: the 3rd user, female, 168cm, 28 years old, can set her parameters as follow:


NOTE:

- During parameter setting, if idle for 6 seconds, the parameters already set will be memorized automatically and the LCD shows "0", meaning that body fat measuring will be started.
- If the LCD shows "0" during the parameter setting, press SET for 3 seconds to return to the setting mode, then you can set the user code, gender, height and age

(2) To Start Measuring (Make sure your feet are touching with the electrodes on the scale platform, or your body fat can not be measured.)

To ensure accuracy, please remove your shoes & socks before stepping on.
 Example: the 3rd user, female, 168cm, 28 years old can measure as follow:

WARNING INDICATIONS

Low Battery Indication
 The battery power is running low, please replace with a new battery.

Over-load Indication
 The weighing subject on the platform exceeds the maximum capability of scale. Please step off to avoid damage.

Low Fat% Indication
 The Fat% is too low. Please lead a more nutritional diet and take good care of yourself.

High Fat% Indication
 The Fat% is too high. Please watch your diet and do more exercise.

Re-test
 Error exists, please re-test to get the correct result.

ABOUT FAT

1. About Fat%

- Fat% is an index % of fat content in human body.
- The fat content in the human body is of a certain level, overabundance of fat will bring bad effects on metabolism, thereby, various illnesses may occur.

2. Working Principles

As electrical signals of a certain frequency transmit through the human body, the impedance of fat is stronger than that of muscle or other human tissues. Through the use of biological impedance analysis, this product will pass a safe frequency of low voltage electrical current through human body, thus, the Fat% will be reflected according to the resistance differences when the electrical signals are transmitting through human body.

3. Fat/Hydration Content Fitness Assessment Chart

After weighing yourself, you can use the following Assessment Chart for reference:

| Age | Fat(F) | Hydration(F) | Fat(M) | Hydration(M) | Status | Key Line Indication |
|-----|------------|--------------|------------|--------------|------------|---------------------|
| <30 | 4.0-16.0% | 66.0-57.8% | 4.0-11.0% | 66.0-61.2% | Thin | |
| | 16.1-20.5% | 57.7-54.7% | 11.1-15.5% | 61.1-58.1% | Slim | |
| | 20.6-25.0% | 54.6-51.6% | 15.6-20.0% | 58.0-55.0% | Normal | |
| | 25.1-30.5% | 51.5-47.8% | 20.1-24.5% | 54.9-51.9% | Heavy | |
| | 30.6-60.0% | 47.7-27.5% | 24.6-60.0% | 51.8-27.5% | Overweight | |
| >30 | 4.0-20.0% | 66.0-55.0% | 4.0-15.0% | 66.0-58.4% | Thin | |
| | 20.1-25.0% | 54.9-51.6% | 15.1-19.5% | 58.3-55.3% | Slim | |
| | 25.1-30.0% | 51.5-48.1% | 19.6-24.0% | 55.2-52.3% | Normal | |
| | 30.1-35.0% | 48.0-44.7% | 24.1-28.5% | 52.2-49.2% | Heavy | |
| | 35.1-60.0% | 44.6-27.5% | 28.6-60.0% | 49.1-27.5% | Overweight | |

The above information is only for reference

ABOUT CALORIE

1. Calorie & Energy Consumption

Calorie: Energy Unit

The data reading indicates the energy you need to sustain the daily metabolism and ordinary activities required for your weight, height, age and sex.

2. Calorie & Losing Weight

Controlling your calorie intake combined with proper exercise is an effective way to lose weight. This is a very simple theory. As the daily energy intake can not support the energy consumption of the human body, the saccharide and fat stored inside will decompose and provide energy to the human body. This will leads to weight loss.

LCD has no response

First check that the batteries are still live, if not, replace them if necessary.

Repairs may only be carried out by our authorized customer service or dealers.

ADVISE FOR USE & CARE

1. DO NOT use the scale in humid, overheat or overcold environment.(Temperature range form 5°C to 35°C)
2. DO NOT use the scale on soft mat,lawn,slippery floor or uneven floor.
3. DO NOT drop,shock or strike the scale as it is a precise instrument.
4. Keep the scale clean.DO NOT use acid cleaning agent.
5. DO NOT use the scale on a place with intense shake. Keep standing firmly on scale when in use.
6. To ensure accuracy, please test at the same time of the sequential days.

SPECIFICATIONS

1. Model: NOPS972
2. 74 x 35,6 mm LCD
3. Low battery / over load indication
4. Capacity: 150 kg / 330 lb / 24st
5. Graduation: 100g / 0.2 lb
6. Power: 1x3V lithium cell(CR2032)
7. Uses a high precision strain gauge sensor
8. Multifunction: tests body fat %, body hydration %,body muscle%, body bone% and suggests calorie intake requirements.
9. Built-in memory for 12 different users
10. Optional normal weighing mode (bypassing parameter setting): Auto step on weighing and consecutive add-on weighings
- 11.Low battery and overload indication

This symbol on the product or in the instructions means that your electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There are separate collection systems for recycling in the EU.

Notes for battery disposal

The batteries used with this product contain chemicals that are harmful to the environment. To preserve our environment, dispose of used batteries according to your local laws or regulations. Do not dispose of batteries with normal household waste.

For more information, please contact the local authority or your retailer where you purchased the product.

Imported and Exclusively marketed by:
Power International AS,
PO Box 523, N-1471 Lørenskog, Norway

Power Norge:
<https://www.power.no/kundeservice/> T: 21 00 40 00
Expert Danmark:
<https://www.expert.dk/kundeservice/> T: 70 70 17 07
Power Danmark:
<https://www.power.dk/kundeservice/> T: 70 33 80 80
Punkt 1 Danmark:
<https://www.punkt1.dk/kundeservice/> T: 70 70 17 07
Expert Finland:
<https://www.expert.fi/tuki/asiakaspalvelu/> T: 020 7100 670
Power Finland:
<https://www.power.fi/tuki/asiakaspalvelu/> T: 0305 0305
Power Sverige:
<https://www.power.se/kundsservice/> T: 08 517 66 000

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