MANO HAND MASSAGER

POMANO
User Manual

GB

Thank you for purchasing your new POINT MANO HAND MASSAGER. These operating instructions will help you use it properly and safely. We recommend that you spend some time reading this instruction manual in order that you fully understand all the operational features it offers. Read all the safety instructions carefully before use and keep this instruction manual for future reference.

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


IMPORTANT INSTRUCTIONS

READ AND SAVE THESE IMPORTANT SAFETY INSTRUCTIONS

1.  Warning: This appliance is not a toy. Keep away from infants, babies and children.
2. This appliance can be used by children aged 8 years and older and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they are supervised or have been instructed on the safe use of the appliance and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be carried out by children without supervision.
3. Children under 8 years old are prohibited from using this product.
4. Do not use this appliance for more than 30 minutes.
5. Do not fall asleep while using this appliance.
6. Pregnant women, infants, children, individuals with swollen or inflamed areas of the body, people suffering from severe muscle or joint pain, and anyone using a pacemaker are not recommended to use this device.
7. Do not use this device if you have taken painkillers, sedatives, or consumed alcohol.
8. Close supervision is required when disabled or ill individuals operate this appliance.


9. Massage should be comfortable and pleasant, if you feel pain or discomfort, discontinue use and consult your doctor.
10. Do not place the appliance near water, under direct sunlight or in excessive heat.
11. Do not use this product when lying flat or reclining.
12. Do not use this product while performing precise tasks or driving.
13. Store this product out of the reach of children.
14. Do not scratch or puncture this product with sharp objects.
15. Do not disassemble the product's internal components, as this may cause damage.
16. This is a household electronic healthcare product, not a medical device, and is not intended for medical diagnosis or treatment.
17. Please use the included charging cable and a safety-compliant power adapter to ensure electrical safety.
18. Excessive use may cause discomfort or skin irritation.

PRODUCT OVERVIEW

USB-C Charging Cable

Safety Air Valve Knob

OPERATING INSTRUCTIONS

⏻ Power Button

Press the power button for 1.5 seconds to turn on the device; the indicator light will stay solid blue. Press the power button for 1.5 seconds to turn the device off.

Ⓜ Mode Button

After turning on the device, the white Ⓜ icon will light up. Briefly press this button to cycle through modes in the order 1 → 2 → 3 → 4 → 5 → 6 → 1, and the corresponding mode indicator will display P1 to P6 in turn.

Pause button

Short press the button to activate the pause function, and the indicator will turn green. Short press again to deactivate the pause function; the device will resume operation, and the indicator light will return to white.

Level button

Short press the button to cycle through six massage intensity levels, and the corresponding intensity level will be displayed.

Timer button

Short press the button to cycle through 10/15/20 minutes.

Heat button

Heat function is disabled by default. Short press the button to activate low heating level, and the indicator will turn light red; short press again to switch to high heating level, and the indicator will turn dark red; short press a third time to turn off the heat function.

Knead button


Short press the button to activate the kneading function, and the indicator light will turn green. Short press again to deactivate the function, and the light will return to white.

Safety air valve knob

If a malfunction occurs and your hand cannot be removed, please rotate the knob to the ON position. If the compression function fails to activate after the device is powered on, please rotate the knob to the OFF position.

Charging

- Fully charge the product for approximately 4.5 hours before the first use.
- Connect the USB-C cable to the product's charging port, then connect the other end to a certified 5V/1A adapter, and plug the adapter into an AC power outlet.

- While charging, the four battery bars will light up from current battery level to four in a circular sequence, indicating that it is charging. When fully charged, all four battery bars will remain lit.
- Low battery warning: The battery indicator will flash in red.

Note:

1. While charging, the product cannot be turned on and used.
2. It is normal for the device and the power adapter to heat up slightly during charging.

Operation

- Visually inspect the product for any shipping damage before use.
- Place your hand inside the massager, and adjust it to the most comfortable position.
- Turn on the massager, then adjust it to your desired mode and intensity level.
- Turn off the massager after use; note that the device will automatically turn off when the timer runs out.

Care and Maintenance

1. Keep the product clean, and keep it away from places with high humidity or high temperature, as well as environments with volatile gases.
2. Protect it from moisture, rodent damage and dirt; do not place any heavy items on top of the product.
3. If the product malfunctions, please seek repair from qualified professionals. Non-professionals must not attempt to disassemble it.
4. Stop using the product immediately if you notice any of the following abnormalities: a burnt smell, intermittent power outages, overheating of the power cord, damaged cord insulation, or any other abnormal conditions.

Troubleshooting

This section summarizes the most common problems you could encounter with the product. If you are unable to solve the problem with the information below, contact the service center immediately.

Problem	Possible cause	Solution
Can not power on	Battery power is insufficient	Charge it, then try again
The product does not work	Low battery	Charge it, then try again
Stops during use	The Timer default working time is up	Long press the power button to restart
	Low battery	Charge it then use it again
Weak massage	Low battery	Fully charge it
No heat compress	The hot compress function is not turned on	Short press the heating button to turn on the hot compress function
The charging indicator light is not on	The USB-C socket has poor contact /poor power contact	Pull and plug the Type-C charging cable / power cord

Specification:

Model: POMANO

Input: DC 5V/1A

Rated power: 5W

Battery: 4000mAh lithium battery

Power consumption: < 0.5 W (Standby)

The device enters standby mode when it is fully charged and the adapter is still connected.

Battery Information

Battery Type: Rechargeable Li-ion Battery Pack

Model Identification: CY 18650-N 4000mAh 1S2P

Battery Category: Portable Battery

Weight: Approx. 95g

Nominal Voltage: 3.7V

Rated Capacity: 4000 mAh / 14.8 Wh

Chemical System: Lithium Nickel Manganese Cobalt Oxide (NMC), Graphite

Hazardous Substances (Hg, Cd, Pb): 0% Mercury, Cadmium, or Lead

Other Hazardous Substances: None

Critical Raw Materials (>0.1%): Lithium, Cobalt, Nickel, Manganese, Graphite, Copper, Aluminium

Extinguishing Agents: Water spray, CO₂, or dry sand

Manufacturer: Xiamen Invcy Electronics Co.,Ltd

Address: 1F/3F-5F, BLDC 5, No.943-5, Tonglong 2nd Rd, Torch Hi-tech Industrial Zone (Xiang'an) Xiamen 361000 Fujian P.R. China.

Origin: Made in China

Notes for battery disposal

The batteries used with this product contain chemicals that are harmful to the environment. To preserve our environment, dispose of used batteries according to your local laws or regulations.

Do not dispose of batteries with normal household waste. For more information, please contact the local authority or your retailer where you purchased the product.

This symbol on the product or in the instructions means that you electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There is separate collection systems for recycling in the EU. For more information please contact the local authority or the retailer where you purchased the product.

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Tillverkare: Xiamen Invcy Electronics Co.,Ltd