EIRA EYE MASSAGER

POEIRA
User Manual

GB

Thank you for purchasing your new POINT EIRA eye massager. These operating instructions will help you use it properly and safely. We recommend that you spend some time reading this instruction manual in order that you fully understand all the operational features it offers. Read all the safety instructions carefully before use and keep this instruction manual for future reference.

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SAFETY INSTRUCTIONS

To ensure safe operation and prevent injury, please read and follow all safety instructions.

Medical and health restrictions

1. Remove contact lenses before using this device.
2. Do not use this device if you are recovering from eye surgery.
3. Do not use this device if you are receiving medical treatment for ocular discomfort.
4. Do not use this device if you suffer from hypertension, high intraocular pressure, arteriosclerosis, or cerebral thrombosis.
5. Do not use this device if you have a persistent high fever, malignant tumors, or osteoporosis.
6. Do not use this device if you have an infectious disease.
7. We do not recommend this device for pregnant women, infants, children, or people with swollen or inflamed body areas, or severe muscle or joint pain.
8. We do not recommend this device for individuals who use a pacemaker.
9. Do not use this device if you have taken painkillers or sedatives, or if you have consumed alcohol.
10. Massage should be comfortable and pleasant. If you feel pain or discomfort, stop using the device immediately and consult your doctor.

11. This device is a household electronic healthcare product. It is not a medical device and is not intended for medical diagnosis or treatment.

Supervision and safety of children

1. Children under 8 years old must not use this device. Keep the device and its packaging away from infants, babies, and children.
2. Children aged 8 years and older, and persons with reduced physical, sensory, or mental capabilities or a lack of experience and knowledge, can use this device under supervision. They must receive instruction on how to use the device safely and understand the hazards involved.
3. Do not allow children to play with the device.
4. Children must not clean or maintain the device without supervision.
5. Store this device out of the reach of children.

Operation and usage guidelines

1. Do not use this device for more than 20 minutes at a time.
2. Do not fall asleep while using this device.
3. Do not use this device while performing precise tasks or driving.
4. Excessive use can cause discomfort or skin irritation.
5. Close supervision is required when disabled or ill individuals operate this device.


Power, charging, and maintenance

1. Do not place this device near water, under direct sunlight, or in excessive heat.
2. Do not scratch or puncture this device with sharp objects.
3. Do not disassemble the internal components of this device, as this can cause damage.
4. Use only the included charging cable and a safety-compliant power adapter to ensure electrical safety.

Product features

1. Voice broadcast.
2. Four built-in music tracks.
3. Bluetooth function for music playback.
4. Next-song adjustment function.
5. Speaker mute option.
6. Vibration function.
7. Two heating levels: 34 °C and 38 °C.
8. Cooling sensation level: 22 °C.
9. Two-button operation with seven modes.
10. Low battery reminder.
11. USB-C charging port.

PRODUCT OVERVIEW

GB

① Power Button:

- **Long press:** Turn the device on or off.
- **Short press:** Cycle through four active modes (Automode, Relax mode, Relieve mode, and Energy mode).
- **Double press:** Switch to the high heat level. The heating indicator turns blue.


- **Default on power-up:** The device starts at the low heat level by default, and the heating indicator turns white.
- ② **Sleep Button:**
 - **Short press:** Cycle through three sleep modes (Soothing mode, Sleep mode, and Cold mode).
 - **Long press:** Mute or unmute the built-in speaker.
 - **Double press:** Play the next music track.
- ③ **Vibration indicator.**
- ④ **Speaker mute and unmute indicator.**
- ⑤ **Bluetooth indicator.**
- ⑥ **Air compression indicator.**
- ⑦ **Heating indicator.**
- ⑧ **Cold indicator.**
- ⑨ **Battery charging indicator.**
 - **Low power (under 15%):** The red LED flashes.
 - **Charging:** The blue LED flashes.
 - **Fully charged:** The blue LED remains on.
- ⑩ **Adjustable strap.**
- ⑪ **Compressing airbag.**
- ⑫ **Hot and cold nodes.**
- ⑬ **USB-C charging port .**

Mode functions

Mode	Air compression	Heat	Cold	Vibration	Speaker and music	Bluetooth
Auto mode	Yes	Yes	Yes	Yes	Yes	Yes
Relax mode	Yes	Yes		Yes	Yes	Yes
Relieve mode	Yes		Yes	Yes	Yes	Yes
Energy mode	Yes				Yes	Yes
Soothing mode		Yes	Yes		Yes	Yes
Sleep mode		Yes			Yes	Yes
Cold mode			Yes		Yes	Yes

Charging

- Fully charge the device for approximately 3.5 to 4 hours before the first use.
- Connect the included USB-C cable to the USB-C charging port of the device. Connect the other end of the cable to a certified 5 V, 1 A power adapter, and plug the adapter into a power outlet.

- The battery charging indicator flashes blue during charging and remains blue when the battery is fully charged.
- Unplug the power adapter and disconnect the charging cable from the device once charging is complete.

- Note that you cannot turn on or use the device while it is charging.
- Note that it is normal for the device and the power adapter to warm up slightly during charging.

Operation

1. Inspect the device for any visible damage before use. Do not use the device if it is damaged.
2. Ensure your face and eye area are clean and dry. Remove any contact lenses before use.
3. Adjust the strap so the device fits snugly and comfortably around your eyes. Make sure there is no excessive pressure on your eyes.
4. Long press the power button to turn on the device, then choose your desired mode.
5. Turn off the device after use by long pressing the power button. The device will also turn off automatically when the timer runs out.

Timer function

To prevent excessive massage, the device features an automatic timer. The device will turn off automatically after 10 minutes of continuous operation.

Bluetooth connection

1. Enable Bluetooth on your mobile device.
2. Search for the eye massager under the name Eira in your Bluetooth settings and pair the devices.

3. Once connected, you will hear a voice prompt saying "Bluetooth connected" and the Bluetooth indicator will light up.
4. You can now play your own music or audio through the built-in speakers of the device.

Care and Maintenance

1. Gently wipe the device with a soft, lint-free cloth or a cotton swab to remove skin oils, skincare products, or sweat.
2. Keep the device clean and dry. Store it away from high humidity, high temperatures, and volatile gases.
3. Protect the device from moisture, dirt, and pests. Do not place heavy items on top of the device.
4. Do not attempt to disassemble or repair the device yourself. Only qualified professionals should service the device.
5. Stop using the device immediately if you notice any problems, such as a burnt smell, overheating, intermittent power loss, or other unusual behavior.

Troubleshooting

This section summarizes the most common problems you could encounter with the device. If you cannot solve the problem using the solutions below, contact customer service immediately.

Problem	Possible cause	Solution
The device does not turn on	The battery is low	Charge the battery and try again
The device does not work	The battery is low	Charge the battery and try again
The device stops working during use	The ten-minute timer has ended	Long press the power button to turn on the device again
The device stops working during use	The battery is low	Charge the battery and try again
The massage feels weak	The battery is low	Fully charge the battery
The air compression does not work	The selected mode does not include air compression	Switch to a mode that includes air compression
The cold function does not work	The selected mode does not include the cold function, or the device has been used for too long	Switch to a mode that includes the cold function, or turn off the device for 30 minutes before using it again
The heat function does not work	The selected mode does not include the heat function	Switch to a mode that includes the heat function
The battery charging indicator does not light up	The USB-C charging cable or the power adapter is poorly connected	Disconnect and reconnect the USB-C charging cable and the power adapter

Specifications:

Model: POEIRA

Input: DC 5V 1A

Rated power: 5W

Power consumption: < 0.5 W (Standby)

The device enters standby mode when it is fully charged and the adapter is still connected.

Battery Information

Battery Type: Rechargeable Li-ion Battery Pack

Model Identification: 104050

Battery Category: Portable Battery

Weight: Approx. 40.8 g

Nominal Voltage: 3.7 V

Rated Capacity: 2500 mAh / 9.25 Wh

Chemical System: Lithium Nickel Manganese Cobalt Oxide (NMC), Graphite

Hazardous Substances (Hg, Cd, Pb): 0% Mercury, Cadmium, or Lead

Other Hazardous Substances: None

Critical Raw Materials (>0.1%): Lithium, Cobalt, Nickel, Manganese, Graphite, Copper, Aluminium

Extinguishing Agents: Water spray, CO₂, or dry sand

Manufacturer: HDY ELECTRONIC TECHNOLOGY CO., LTD.

Address: Floor 6, Building 1, Jiangnan Industrial Building, No. 88, Jiangnan South, Xiangang Community, Chang'an Town, Dongguan, China

Notes for battery disposal

The batteries used with this product contain chemicals that are harmful to the environment. To preserve our environment, dispose of used batteries according to your local laws or regulations.

Do not dispose of batteries with normal household waste. For more information, please contact the local authority or your retailer where you purchased the product.

This symbol on the product or in the instructions means that you electrical and electronic equipment should be disposed at the end of its life separately from your household waste. There is separate collection systems for recycling in the EU. For more information please contact the local authority or the retailer where you purchased the product.

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Tillverkare: HDY ELECTRONIC TECHNOLOGY CO., LTD.